

IQ10 Accessories

TROLLEYS

KEN HYGIENE SYSTEMS®

iQ¹⁰



IQ10 Cart Washer

Product

The IQ10 is a Cart Washer designed to meet diverse cleaning and disinfection needs. It is engineered to deliver low consumption of electricity, water, and detergents. Its internal water recirculation system collects the disinfection water and reuses it, reducing overall water usage.

With a short process time, the IQ10 offers a time-efficient washing solution, making it particularly useful for environments where maintaining a high throughput is important.

Especially suitable for cleaning and disinfecting larger equipment due to its spacious design and adaptable trolleys. It is ideal for washing surgical instruments, utensils, vessels, and containers.

Additionally, it is a perfect solution for care facilities to wash wheelchairs, toilet chairs, walking aids, bowls, and similar equipment.

The IQ10 complies with the latest regulatory standards, including EN/ISO 15883-1/2/6. Its compliance allows it to serve as a backup for standard Washer Disinfectors in a CSSD.



IQ10 Cart Washer

Four sizes

IQ10: Standard version

External dimensions: H = 2790 mm, W = 2540 mm, L = 2770 mm

Washing chamber dimensions: H = 2080 mm, W = 1250 mm, L = 2350 mm

IQ10-W: Wide version

External dimensions: H = 2790 mm, W = 2940 mm, L = 2770 mm

Washing chamber dimensions: H = 2080 mm, W = 1450 mm, L = 2350 mm

IQ10-L: Long version

External dimensions: H = 2790 mm, W = 2540 mm, L = 3520 mm

Washing chamber dimensions: H = 2080 mm, W = 1250 mm, L = 3100 mm

IQ10-LW: Long and wide version

External dimensions: H = 2790 mm, W = 2940 mm, L = 3520 mm

Washing chamber dimensions: H = 2080 mm, W = 1450 mm, L = 3100 mm



Container Trolleys

Trolley for containers and lids, which can be loaded from both sides. Each shelf is designed to hold multiple containers. The shelves are easily removed and can be moved to other positions to accommodate a large range of different container.

The trolleys can be delivered in two lengths:

- **Version A. ref. 2992114.**
Dimensions: Length: 1061 mm. Height: 1900 mm. Depth: 1020 mm.
Weight: 75 kg.
Capacity: 12 containers and room for 12 lids on the top shelf.
Fits in IQ10 and IQ10-W.
- **Version B. ref. 2990966.**
Dimensions: Length: 1436 mm. Height: 1900 mm. Depth: 1020 mm.
Weight: 87 kg.
Capacity: 16 containers and room for 16 lids on the top shelf.
Fits in IQ10-L and IQ10-LW.

Optionally, a shelf for lids can be chosen on the top. One lid holder can contain 4 lids.

- **Holder for lids. ref. 2998509.**
Dimensions: 615 mm. x 279 mm.

Depending on the length of the trolley and the machine, it is possible to place two trolleys in the washing chamber.



Instrument Trolleys

4-Level instrument washing trolley, loadable from both sides. With rotating spray arms positioned between each shelves to ensure direct spray on the material placed on the shelves. This improves the result and efficacy of the cleaning process.

The shelves and spray arms are easy to remove.

The trolleys can be delivered in 2 lengths:

- **Version A.** ref. 2992101.
Dimensions: Length: 1061 mm. Height: 1900 mm. Depth: 1020 mm.
Weight: 85 kg.
Capacity: 16 DIN baskets.
Fits in IQ10 and IQ10-W.
- **Version B.** ref. 2000965.
Dimensions: Length: 1436 mm. Height: 1900 mm. Depth: 1020 mm.
Weight: 95 kg.
Capacity: 32 DIN baskets.
Fits in IQ10-L and IQ10-LW.

Depending on the length of the trolley and the machine, it is possible to place two trolleys in the washing chamber.



Trolleys for OP-Shoes

Trolley suitable for OP shoes, and other shoes. It can be loaded from both sides. The shelves can easily be removed and moved to other positions.

The trolleys can be delivered in two lengths:

- **Version A.** ref. 2992121.
Dimensions: Length: 1061 mm. Height: 1900 mm. Depth: 1020 mm.
Weight: 77 kg.
Capacity: 70 shoes.
Fits in IQ10 and IQ10-W.
- **Version B.** ref. 2990967.
Dimensions: Length: 1436 mm. Height: 1900 mm. Depth: 1020 mm.
Weight: 91 kg.
Capacity: 100 shoes.
Fits in IQ10-L and IQ10-LW.

Depending on the length of the trolley and the machine, it is possible to place two trolleys in the washing chamber.



Special Trolley

Special trolley which can be change to suit the items needed to be washed. The shelves can easily be removed and moved to other positions.

This trolley is suitable for washing wheel chairs, and other similar devices and aids.

Only one version. ref. 2992085

Dimensions: Length: 1588 mm. Height: 1900 mm. Depth: 1020 mm.
Weight: 77 kg.

Depending on the length of the trolley and the machine, it is possible to place two trolleys in the washing chamber.



Malone Hospital Services Pty Ltd.
24 Tathra Street West Gosford, 2250, NSW
Phone: +61 (02) 4323 9731
sales@malones.com.au
www.malones.com.au



Disclaimer

Do not use this technical data sheet for installation of equipment. Subject to change without notice. MEJ · 01.25

KEN HYGIENE SYSTEMS A/S · Bøgebjergvej 60 · DK-5672 Broby · Phone +45 6263 1091 · ken@ken.dk · www.ken.dk